

# Coverage® Spray TB Plus

## 21% Isopropyl Alcohol and Quaternary Ammonia blend Broad Spectrum Healthcare Surface Disinfectant and Sanitizer

### PRODUCT DESCRIPTION

Virucidal, Antimicrobial, and Fungicidal surface disinfectant for use in healthcare settings

### WHERE TO USE

Coverage® Spray TB Plus can be used in healthcare facilities.

FORMAT: LIQUID, Ready-To-Use



### DIRECTIONS FOR USE:

It is a violation of Federal law to use Coverage® Spray TB Plus in a manner inconsistent with its labeling.

Pre-clean all surfaces to be treated with Coverage® Spray TB Plus. Apply this undiluted product by spray to wet surfaces thoroughly.

**For surface cleaning:** To clean surfaces, saturate surfaces with Coverage® Spray TB Plus by spray then wipe surface clean. Repeat, if necessary, until surfaces are visibly clean. Discard paper towels or launder cloth towels before reusing.

**For cleaning walls and floors:** In a clean bucket, dilute 2 to 4 oz. of Coverage® Spray TB Plus per gallon of warm, clean, potable water. Thoroughly damp-mop using a clean mop, cloth towel, or sponge. Allow to air dry.

**For use as a pre-cleaning Instrument spray:** Place item into a suitable container. Thoroughly spray item with Coverage® Spray TB Plus until drenched, thoroughly rinse. Follow with an appropriate terminal cleaning and sterilization/disinfection process.

**For hard, non-porous, non-food contact surface sanitization:** Apply Coverage® Spray TB Plus to pre-cleaned hard, non-porous surfaces by spray to wet surfaces thoroughly. Allow hard, non-porous surfaces to remain visibly wet for 10 seconds at room temperature then wipe surfaces using a clean mop, paper or cloth towel, or rinse surface using potable water or allow to air dry. Discard or launder spent materials before reuse.

**Note:** Food contact surfaces must be rinsed with potable water after sanitization and before reuse.

**For use as a non-critical medical item disinfection immersion solution:** Thoroughly pre-clean items being sure to remove all gross filth and heavy soil, then rinse and dry. Allowing for complete submersion of item(s), place pre-cleaned item(s) into an appropriate container that has been filled with Coverage® Spray TB Plus. Flush lumens of hollow devices with Coverage® Spray TB Plus. Allow item(s) to remain submerged for 2 minutes. Remove, rinse with clean, potable water and dry. Change solution daily or when it becomes visibly soiled or discolored.

**For use as a hard, non-porous non-food contact surface disinfectant/virucide/fungicide:** Apply Coverage® Spray TB Plus to pre-cleaned surfaces by spray, to wet surfaces thoroughly. Allow to remain wet for 2 minutes at room temperature, then wipe dry using a mop, clean paper or cloth towel, or rinse surface using potable water and wipe surface dry or allow to air dry. Discard or launder spent materials before reuse.

Coverage® Spray TB Plus is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. Coverage® Spray TB Plus may be used to pre-clean critical and semi-critical medical devices prior to sterilization or high-level disinfection.

# Coverage® Spray TB Plus

## Broad Spectrum Healthcare Surface Disinfectant and Sanitizer

### KILL CLAIMS:

#### Hospital and General Disinfection or Bactericidal:

Coverage® Spray TB Plus kills the following bacteria on pre-cleaned hard, non-porous surfaces in 2 minutes at room temperature when used as directed:

Organism	Kill time
Acinetobacter baumannii ATCC 19606	2 minutes
Escherichia coli -with (extended beta-lactamase resistance ESBL) (E. coli) ATCC BAA-196	
Enterobacter aerogenes ATCC 13048	
Carbapenem-resistant Klebsiella pneumoniae (CRKP) ATCC BAA-1705	
Listeria monocytogenes ATCC 19117	
Methicillin-Resistant Staphylococcus aureus (MRSA) ATCC 33591	
Mycobacterium bovis BCG (tuberculosis) (TB) Organon Teknika Corp	
Pseudomonas aeruginosa ATCC 15442	
Serratia marcescens ATCC 14756	
Salmonella enterica ATCC 10708	
Staphylococcus aureus ATCC 6538	
Vancomycin-Resistant Enterococcus faecalis (VRE) ATCC 51559	

#### Soft, Non-Food Contact Surface Sanitization

Coverage® Spray TB Plus kills the following organisms in 10 seconds on pre-cleaned soft surfaces at room temperature when used as directed:

Organism	Kill time
Enterobacter aerogenes ATCC 13048	10 seconds
Staphylococcus aureus ATCC 6538	

#### Fungicidal

Coverage® Spray TB Plus kills the following fungi in 2 minutes on pre-cleaned hard, non-porous surfaces at room temperature when used as directed:

Organism	Kill time
Trichophyton interdigitale (formerly known as Trichophyton mentagrophytes) (Ringworm) (Athletes Foot) ATCC 9533*	2 minutes

#### Hard, Non-Porous Non-Food Contact Surface Sanitization

Coverage® Spray TB Plus kills the following organisms in 10 seconds on pre-cleaned hard, non-porous non-food contact surfaces at room temperature when used as directed:

Organism	Kill time
Enterobacter aerogenes ATCC 13048	10 seconds
Staphylococcus aureus ATCC 6538	

#### Virucidal\*

Coverage® Spray TB Plus kills the following viruses\* on pre-cleaned hard, non-porous surfaces in 2 minutes at room temperature when used as directed:

Organism	Kill time
Avian Influenza A Reassortant Virus (H3N2) ATCC VR-2072	2 minutes
Influenza A Virus (H1N1) Strain A/Swine/Iowa/15/30 (formerly known as "swine flu") ATCC VR-333	
Rotavirus Strain WA	
Hepatitis B Virus* (HBV) Strain LeGarth	
Hepatitis C Virus* (HCV)	
Herpes Simplex Virus (HSV-2) ATCC VR-734	
Human Immunodeficiency Virus* (HIV-1) (associated with AIDS) HTLV-111 B	
Influenza A Virus Strain A2/Hong Kong ATCC VR-544	
Rhinovirus Type 14 RV-14 ATCC VR-284	
SARS-Related Coronavirus 2 SARS-CoV-2 USA-WA1/2020	

### \*HIV-1, HBV AND HCV PRECAUTIONS

Coverage® Spray TB Plus kills HIV-1, HBV and HCV on pre-cleaned hard, non-porous surfaces/objects previously soiled with blood/body fluids in health care or other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1)(associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV).

#### Special instructions for cleaning against HIV-1, HBV and HCV on surfaces/objects soiled with blood/body fluids:

**Personal Protection:** Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks and eye coverings.

**Cleaning Procedure:** Blood and body fluids must be thoroughly cleaned from surfaces/objects before disinfection.

**Contact Time:** Allow hard, non-porous surfaces to remain wet for 2 minutes to kill HIV-1, HBV and HCV.

**Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids must be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

# Coverage<sup>®</sup> Spray TB Plus

## Broad Spectrum Healthcare Surface Disinfectant and Sanitizer

### STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

**Pesticide Storage:** Store in original container in a cool, well-ventilated (preferably locked) area inaccessible to children and pets. **Pesticide Disposal:** To avoid waste, use all material in this container by application according to label directions. If wastes cannot be avoided, offer remaining product to a waste disposal facility or pesticide disposal program. **Container Handling:** Non-refillable container. Do not reuse or refill this container. Discard container and put in trash or offer for recycling if available.

### ACTIVE INGREDIENTS

N-ALKYL (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> )	
DIMETHYL BENZYL AMMONIUM CHLORIDE .....	0.154%
N-ALKYL (68% C <sub>12</sub> , 32% C <sub>14</sub> )	
DIMETHYL ETHYLBENZYL AMMONIUM CHLORIDE.....	0.154%
ISOPROPANOL .....	21.000%
<b>OTHER INGREDIENTS</b> .....	<u>78.692%</u>
<b>TOTAL</b> .....	100.000%

### PRECAUTIONARY STATEMENTS

#### Hazards to Humans and Domestic Animals

**CAUTION:** Harmful if inhaled or absorbed through skin. Causes moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wear appropriate personal protective equipment based on safe work practices. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

### PHYSICAL OR CHEMICAL HAZARDS

Do not use or store near heat or open flame.

### FIRST AID

If inhaled: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact the National Pesticides Information Center (NPIC) at 1-800-858-7378 for non-emergency information concerning Coverage<sup>®</sup> Spray TB Plus.



### PACK SIZES

Stock Code	Size	Case Qty
1629B4	24 oz Spray Bottle	12
162908	1 Gallon Bottle	4

SC Johnson Professional USA, Inc.  
2815 Coliseum Center Dr.,  
Suite 600, Charlotte, NC 28217  
United States  
Telephone: +1 866-783-0422

[www.scjp.com/healthcare](http://www.scjp.com/healthcare)

For Hospital and Professional Use Only. This information and all further technical advice is based upon our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. For external use only. Keep out of reach of children. If swallowed, get medical help or contact a Poison Control Center right away.

