

Cutan Optidose Antiseptic Handrub

PRODUCT EDUCATION SHEET



DESCRIPTION: An antiseptic handrub that combines unique¹ alcohol foam with added conditioners to help hospitals keep patients safe and improve compliance with just **ONE PUSH**.

USAGE: For use in healthcare environments where high levels of hand hygiene are vital. Can be used for surgical hand disinfection.

METHODS OF USE: Gently press once fully to apply product in to cupped hands, covering all surfaces. Rub into hands using the WHO recommended hand hygiene technique for at least 30 seconds. A Cutan Optidose handrub technique poster is available. Repeat as needed throughout the day.

SHELF LIFE:

This product has a shelf life of at least 48 months from the date of manufacture when stored unopened at room temperature.

Feature	Benefit
One Push	Only one push is required to effectively cover the hand with enough product to achieve the recommended 30 second wet time.
Surgical Scrub	Approved for use as a waterless surgical scrub to the EN12791 standard.
Instant Foam Action	Rapidly and easily rubs across hands to kill most common germs present on the skin.
Dispensed as a foam	Provides user with complete control of the product in the hands, preventing the dripping and splashing associated with using liquid/gel alcohol sanitisers.
Convenient	Available in a range of pack sizes, including the hygienic ultrasonically sealed 1 litre cartridge for use in Deb Hand Sanitiser dispenser; plus portable pump packs for bed-end systems, personal use and areas of limited space.
Contains moisturisers	Helps prevent skin dryness and leaves the skin feeling soft and refreshed after use.
No Mess	Does not clog pumps or splash and drip on surrounding floor and wall surfaces as associated with liquid and gel sanitisers
No water required	Quickly dries on hands without the need for rinsing with water.

AVAILABLE PACK SIZES

Stock Code	Description	Case Quantity	Dispenser
COP1L	1 Litre Cartridge	6	DIS2160
COP400ML	400mL Pump Bottle	12	Fits all Cutan brackets

¹Australian Patent No. 2004275900
 AUST R 302680